

CENTURION®

# HubGuard®

## CATHETER SECUREMENT



Available individually sterile packaged or as part of complete Dressing Change or Securement Kits

Meets INS Standards for Catheter Stabilization

- Adhesive-backed foam securely anchors catheter at infusion site
- Cushion protects skin and provides comfort from catheter hubs
- Enhances stability to reduce "pistoning" of catheter in vein
- Unlike gauze, HubGuard won't absorb blood
- Non-irritating adhesive is gentle on skin, yet stays in place for 72 hours or longer
- Available in two latex-free materials:
  - For maximum securement and cushioning, 1/16" thick PE foam
  - For maximum breathability, 1/32" thick breathable PVC foam



**CENTURION**  
MEDICAL PRODUCTS


301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843  
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | [www.centurionmp.com](http://www.centurionmp.com)

CENTURION®

# HubGuard®

CATHETER SECUREMENT

## HUBGUARD SIZES AVAILABLE



PE = 1/16" thick Polyethylene

PVC = 1/32" thick Polyvinyl Chloride

\* Release tabs not included in dimensions

Other sizes and styles may be available.

See your Centurion sales representative for additional catheter securement products.

**CENTURION**  
MEDICAL PRODUCTS

301 CATRELL DRIVE | P.O. BOX 170 | HOWELL | MICHIGAN | 48843  
phone 517.546.5400 | toll free 800.248.4058 | fax 517.546.9388 | [www.centurionmp.com](http://www.centurionmp.com)

0623-09