

CENTURION®

CENTRAL LINE INSERTION BUNDLE

Complete kit for
central line insertions



The Centurion Central Line Insertion Bundle gives you everything you need to minimize central line infection risk, and lets you choose those components that best meet your hospital-specific needs.

Featuring:

NEW Edwards™ Pressure Injectable, Antimicrobial Catheters

Centurion SorbaView® SHIELD

BIOPATCH® Protective Disk with CHG

Centurion Sutureless Catheter Securement Devices

Centurion SnagFree® Needle Holders

ChloroPrep®

Latex Free Components and Packaging

Sharps Safety Products

Ultrasound Gel and Probe Covers

Full-body drapes, gowns, masks and other components necessary for optimum protection against CRBSI (Catheter Related Bloodstream Infections)



EDWARDS PRESSURE INJECTABLE CATHETERS

The Vantex Catheter features Oligon Agent which provides broad spectrum antimicrobial protection.

The Vantex has an infusion rate up to 77% faster than competing catheters.

Vantex and Multi-Med Catheters are available with Edwards AMC THROMBOSHIELD, an antimicrobial Heparin coating.



CENTURION CVC SECUREMENT ANCHOR

Centurion CVC Securement Anchor allows for easy application.



CENTURION CENTRAL LINE BUNDLE

The Centurion Central Line Bundle organizes all of your components in an easy to use sequence. Drapes, gowns and caps are conveniently placed in order of use. Pre-filled sterile syringes are packaged sterile-field ready inside the kit. Products are uniquely packaged to your protocol specifications.

SORBAVIEW® SHIELD

SorbaView SHIELD is a one-step catheter securement system that combines the features and benefits of a SorbaView Dressing with SHIELD Technology. SHIELD Technology guards the catheter from aggressive movements and eliminates the need for a secondary securement device. With the addition of SHIELD Technology, dressings stay on and catheters stay put. By keeping the catheter in place, the risk of complications is dramatically reduced.



SORBAVIEW® ULTIMATE IJ DRESSING SYSTEM






Particularly suited to large introducer catheters, the SorbaView Ultimate IJ Dressing System with Tubing Anchor will not fall off your patient. It provides a proven barrier to liquids and bacteria, while allowing excellent breathability. SorbaView is easy to apply, minimizes patient discomfort, and remains intact for extended periods.



SORBAVIEW® ULTIMATE WITH ADHESIVE-FREE ZONE

Our specially designed SorbaView Ultimate Dressing for central lines features an adhesive-free zone that is designed to fit over the BIOPATCH® Protective Disk. This system increases the area of protection for the insertion site and allows for easy removal of the dressing without jeopardizing the integrity of the BIOPATCH.





DRAPES FOR CENTRAL LINE INSERTIONS

In addition to offering a large variety of stock drapes, Centurion can fully customize a drape to suit your specific needs. Many options are available, including:

Code Number: DC4065
 Size: 40" x 65"
 Material: Nonwoven
 Fenestration: 4" with adhesive border
 Perforation: Yes
 Absorbent Zone: 32" x 33"

Code Number: DC65104P
 Size: 65" x 104"
 Material: Nonwoven
 Fenestration: 3¼" with adhesive border
 Perforation: Yes
 Window: 18" x 18" clear plastic

Code Number: DCT1000
 Size: 60" x 90"
 Material: SMS
 Fenestration: 3" x 5" with tearaway to open to 7" x 9"

Code Number: DC6590
 Size: 65" x 90"
 Material: Nonwoven
 Fenestration: 4" with adhesive zone

Code Number: DC72902FP
 Size: 72" x 90"
 Material: SMS and Absorbent/Light
 Fenestration: 2 each, 4" x 5" with adhesive border
 Perforation: 2 each, fenestration to side
 1 each, head opening


Contact your representative for more options.

Material Key:



PROTECTIVE GEAR

- Full-body drapes
- Surgical gowns – non-reinforced, fabric reinforced and impervious – with multiple protection levels to choose from
- Additional gowns available for support clinicians and extra personnel present during procedure
- Gloves
- Caps and bouffants
- Face masks, with or without splash shields



ULTRASOUND GEL AND PROBE COVERS



STERILE LABELS

Sterile labels can be custom made to fit your exact procedural needs to enable compliance with medication labeling requirements.






SALINE, HEPARIN AND LIDOCAINE SOLUTIONS

A large selection is available.



PRE-FILLED STERILE FIELD READY SYRINGES

Packaged within kits – a Centurion exclusive.



CHLORAPREP APPLICATORS

Chlorhexidine-based applicator is available in 3ml, 10ml, tinted and non-tinted versions.




SNAGFREE® NEEDLE HOLDER

The Centurion SnagFree Needle Holder has a patented design to prevent the suture from snagging on the instrument while tying knots. SnagFree Instruments are precision-designed to provide reusable quality at a disposable price, and are guaranteed to perform to the highest standard every single time.

High-quality scissors and forceps are also available.

HOSPITAL FORMS AND PATIENT/FAMILY EDUCATIONAL MATERIALS

Customized procedural checklists, labels, signs and forms can be built into your trays to ensure proper protocol.



SAFETY COMPONENTS

A large selection of safety scalpels and needles from major manufacturers are available.



VALVES

A large selection of safety valves and closed valve systems are available.



IHI RECOMMENDATIONS

“For the operator placing the central line and for those assisting in the procedure, maximal barrier precautions means strict compliance with hand hygiene and wearing a cap, mask, sterile gown, and gloves. The cap should cover all hair and the mask should cover the nose and mouth tightly. These precautions are the same as for any other surgical procedure that carries a risk of infection.

For the patient, applying maximal barrier precautions means covering the patient from head to toe with a sterile drape, with a small opening for the site of insertion.”

- *Institute for Healthcare Improvement,
100,000 Lives Campaign,
“Prevent Central Line Infections: How-to Guide”*

CDC GUIDELINES

“Maximal sterile barrier precautions (e.g., cap, mask, sterile gown, sterile gloves, and large sterile drape) during the insertion of CVCs substantially reduces the incidence of CRBSI compared with standard precautions...”

- *Centers for Disease Control and Prevention
“Guidelines for the Prevention of Intravascular
Catheter-Related Infections,” CDC Recommendations
and Reports, August 9, 2002 / 51(RR10);1-26*

JCAHO REQUIREMENTS

“Medications and solutions both on and off the sterile field should be labeled even if there is only one medication being used.”

- *Joint Commission on Accreditation of Healthcare Organizations,
2006 Patient Safety Goals*

Centurion, SorbaView and SnagFree are registered trademarks of Centurion Medical Products Corporation. SorbaView® SHIELD and SHIELD Technology are the proprietary property of Centurion Medical Products.

Vantex, Multi-Med and AMC THROMBOSHIELD are trademarks of Edwards Lifesciences Corporation.

BIOPATCH is a registered trademark of ETHICON, INC.

Oligon is a trademark of Implemed, Inc.

ChloroPrep is a registered trademark CareFusion or one of its subsidiaries.

CENTURION
MEDICAL PRODUCTS

100 CENTURION WAY | P.O. BOX 510 | WILLIAMSTON | MICHIGAN | 48895
517.546.5400 | 800.248.4058 | fax 517.546.9388 | www.centurionmp.com