

CENTURION®

# HUBGUARD®

Catheter Securement



## Cushion, Secure and Protect

The Centurion HubGuard is made of an adhesive-backed foam that securely anchors the catheter at the infusion site, while providing a soft barrier between the rigid hub and the patient's skin.

Introducing the new HGDBS – the latest addition to the Centurion HubGuard line. Small enough to be easily used with our small SorbaView® SHIELD (SV254UDT), its double-sided adhesive provides temporary hold of the catheter hub while the dressing is being applied.

## CAREGIVER BENEFITS

- Made with adhesive-backed foam that securely anchors the catheter at infusion site
- Enhances stability to reduce unscheduled I.V. re-starts
















## HOSPITAL BENEFITS

- Enhances stability to reduce "pistoning" of catheter in vein, reducing catheter migration and dislodgement
- Meets INS Standards for Catheter Stabilization

## PATIENT BENEFITS

- Cushion protects skin and provides comfort from rigid catheter hubs
- Available in two latex-free materials:
  - For maximum securement and cushioning, 1/16" thick PE foam
  - For maximum breathability, 1/32" thick breathable PVC foam
- Non-irritating adhesive is gentle on skin, yet stays in place for 72 hours or more

## HUBGUARD SIZES AVAILABLE


<p><b>HGL1</b></p>  <p>3" x 3/4" (PE)</p>	<p><b>HGL3</b></p>  <p>3 3/8" x 3/4" (PVC)</p>	<p><b>HGL4</b></p>  <p>3 3/8" x 3/4" (PE)</p>	<p><b>HGL6</b></p>  <p>1 - 2 3/8" x 5/8"* 1 - 7/8" x 3/8"* (PE)</p>	<p><b>HGDBS</b></p>  <p><b>new</b> 1 - 2 3/8" x 5/8"* 1 - 7/8" x 3/8"* (PE)</p>	<p><b>HG2</b></p>  <p>2 - 1 1/2" x 3/8" (PE)</p>	<p><b>HG10</b></p>  <p>1 - 3 3/8" x 3/4" 3 - 3 3/8" x 1/2" (PE)</p>	<p><b>HGL2</b></p>  <p>2 3/8" x 5/8" (PVC)</p>	<p><b>HG05</b></p>  <p>2 3/8" x 5/8" (PE)</p>	<p><b>HG06</b></p>  <p>1" x 1 1/2" (PE)</p>
<p><b>HG02</b></p>  <p>1 - 3" x 3/4" 1 - 3 3/8" x 3/4" 1 - 1" x 1 1/2" (PE)</p>	<p><b>HG07</b></p>  <p>1 - 3" x 3/4" 1 - 3 3/8" x 3/4" (PE)</p>	<p><b>HG08</b></p>  <p>2 - 3" x 3/4" (PE)</p>	<p><b>HG3</b></p>  <p>1 - 3 3/8" x 3/4" 1 - 1" x 1 1/2" (PE)</p>	<p><b>HG03</b></p>  <p>1 - 3 3/8" x 3/4" 3 - 3 3/8" x 1/2" 1 - 1" x 1 1/2" (PE)</p>					

PE = 1/16" thick Polyethylene

PVC = 1/32" thick Polyvinyl Chloride

\* Release tabs not included in dimensions

Other sizes and styles may be available.  
See your Centurion sales representative for additional catheter securement products.

 Shading denotes double-sided adhesive area

**CENTURION**  
MEDICAL PRODUCTS

100 CENTURION WAY | P.O. BOX 510 | WILLIAMSTON | MICHIGAN | 48895  
517.546.5400 | 800.248.4058 | fax 517.546.9388 | www.centurionmp.com

